


F130

FIXED POOL LIFT

The F130 fixed lift allows, through its flagged structure, to be fixed to perimeter or bearing walls, thus eliminating the need for a basement on the ground. The anchorage must always be carried out on a bearing wall in stone or with a thickness of more than 25 cm by means of dowels or safety counter-plate. The lift can be equipped with a seat + wheels or a sling. The model F130 is particularly suitable for raised therapeutic pools in masonry, whirlpools, mud therapy that are in confined spaces.





Regulation EU 2017/745 Medical Devices – class 1 - DM 2563151/R
Directive 2006/42/EC - Machinery Directive
Directive 2004/108/EC - EMC Electromagnetic Compatibility



Maximum capacity: 140 kg
Maximum travel: 1400 mm (*1400)
Time of lowering - lifting with weight into water: 20 sec
Weight of the lift: 50 kg
Battery: 24 V

Different arms lengths available
Protective cover for electric actuator
Small arch with fabric sling available



TECHNICAL FEATURES

- Seat with removable trolley
- IPX06 spiral cable push-button panel for up-down controls from the water
- Retractable footrest
- Adjustable armrests
- 360° manual rotation
- Removable battery with charger
- Emergency stop button
- Average autonomy: 50 working cycles
- 2-point safety belt

- Galvanizing and powder coating treatment
- Stainless steel frame for in-water parts
- Battery charge level indicator